CHROMATIC COLI/ COLIFORM

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DESCRIPTION

Chromogenic selective medium for *E. coli* and coliforms isolation and enumeration in foods and water. Its particular formulation allows also to carry out the direct indole test to confirm *Escherichia coli* by adding a drop of Kovacs' Reagent (ref.80270) and observing the reagent turning to red.

METHOD PRINCIPLE

Peptone and yeast extract are a source of amino acids and vitamins. The selectivity is due to the presence of the bile salts in the chromogenic mix which inhibit the growth of Gram-positive microorganisms. Beside the chromogenic mix allows the differentiation of *Escherichia coli* from the other enterobacteria. Sodium chloride maintains the osmotic balance of the medium.

TEST PROCEDURE

- Inoculate the surface of the medium with 0.1 ml or 1.0 ml of the specimen to examine.
- Incubate at: 36±1°C for 24 hours for *E. coli*, 44±/0.5°C for 24 hours for faecal coliforms, 30°C for 24-48 hours for total coliforms.

INTERPRETING RESULTS

Interpret results as indicated in table n°1.
Table n°1

Microorganisms	Growth	Typical appearance of the colonies
<i>Escherichia coli</i>	Good	Green
<i>K. pneumoniae/ E. cloacae</i>	Good	Mauve
<i>Proteus mirabilis</i>	Good	Colorless
<i>Pseudomonas aeruginosa</i>	Good	Yellowish
<i>Staphylococcus aureus</i>	Inhibited	
<i>Enterococcus faecalis</i>	Inhibited	

Final identification must be performed by biochemical and/or serological tests.

QUALITY CONTROL

1. Control of general characteristics, label and print
2. Sterility control
7 days at 25 ± 1°C, aerobically
7 days at 36 ± 1°C, aerobically
3. Microbiological control
Plates are inoculated with the following microbial strains: *Escherichia coli* ATCC 25922, *Klebsiella pneumoniae* ATCC 13883, *Enterobacter cloacae* ATCC 23355, *Proteus mirabilis* ATCC 25933, *Pseudomonas aeruginosa* ATCC 27853, *Staphylococcus aureus* ATCC 25923, *Enterococcus faecalis* ATCC 19433.
Inoculum for productivity: 10-100 UFC/ml
Inoculum for specificity: ≤ 10⁴ UFC/ml
Inoculum for selectivity: 10⁴-10⁵ UFC/ml
Incubation conditions: 18-24 h at 36±/1°C.

STORAGE

2-8°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product cannot be classified as hazardous under current legislation nor does it contain harmful substances in concentrations ≥ 1%. It therefore does not require a Safety Data Sheet to be available.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAFIA / BIBLIOGRAPHY

1. Comparison of enumeration of *E. coli* on CHROMagar *E. coli* and MPN methods. Study by S. Weissman, Israel, July 1994.
2. Quantitative determination of *Escherichia coli* in water using CHROMagar *E. coli*. Jose L. Alonso et al. Journal of Microbiological Methods, 25, 1996, p.309-315.

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