# CHROMATIC E.COLI O157

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## BIBLIOGRAFIA / BIBLIOGRAPHY

- L.M. Durso and J.E. Keen USDA, ARS, US Meat Animal Research Center, Clay Center, NE, USA Publication : Journal of Applied Microbiology ISSN 1364-5072: 2007.

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Codice / Code: 11610	Contenuto / Content: 20 piastre / plates
----------------------	--

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## DESCRIPTION

Chromogenic medium for *Escherichia coli* O157 identification from clinical specimens and foods.

## METHOD PRINCIPLE

Peptones and yeast and meat extracts are a source of amino acids, proteins and vitamins. Sodium chloride maintains the osmotic balance of the medium. Sodium Desoxycholate has an inhibitory action. The chromogenic mix allows the identification of the microorganism on the basis of the color and morphology of the colonies. Agar is the solidifying agent.

## TEST PROCEDURE

Inoculate the surface of the medium with the specimen to examine, using a sterile loop.  
Incubate at 36+/-1°C for 18-24 hours. Observe the growth and color of the colonies.

## INTERPRETING RESULTS

Interpret results as indicated in table n°1.  
Table n°1

Microorganisms	Growth	Typical appearance of the colonies
<i>Escherichia coli</i> O157	Good	Mauve or mauve- purple
<i>Escherichia coli</i>	Good	Blue
<i>Klebsiella pneumoniae</i>	Good	Blue
<i>Staphylococcus aureus</i>	Inhibited	
<i>Enterococcus faecalis</i>	Inhibited	

Some Salmonella spp could give false positive results. Final identification must be confirmed by serological tests.

## QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
7 days at 20 ± 2°C  
7 days at 35 ± 2°C
- Microbiological control  
Plates are inoculated with the following microbial strains: *Escherichia coli* O157 H7 ATCC 35150, *Escherichia coli* ATCC 25922, *Klebsiella pneumoniae* ATCC 13883, *Staphylococcus aureus* ATCC 25923, *Enterococcus faecalis* ATCC 19433.  
Incubation conditions: 18-24 h at 35 ± 2°C.

## STORAGE

2-8°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

## WARNING AND PRECAUTIONS

The product cannot be classified as hazardous under current legislation nor does it contain harmful substances in concentrations ≥1%. It therefore does not require a Safety Data Sheet to be available. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

## DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

Liofilchem s.r.l.

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330  
Web site: <http://www.liofilchem.net> E-mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

Rev. 0 /09.02.2011