

CLINET® 101

Chlorhexidine cleaning wipes

Technical notice N° 2160 A (GB)

The finishing touch

HYGIENE

All surfaces and frequently used objects need regular in-depth cleaning. CLINET® 101 disinfecting wipes provide the finishing touch, cleaning and disinfecting at the same time. Give it a wipe, CLINET® takes care of the rest.

EFFICIENCY

The finishing touch is also the touch that finishes bacteria, fungi and viruses, including Herpes, Hepatitis B and HIV 1 (in vitro) virus ; the disinfecting action of CLINET® 101 has been officially tested and certified. CLINET® 101 wipes clean the unseen.

EVERYWHERE

In the laboratory, in the office, at home, CLINET® 101 disinfecting wipes can be used anywhere and everywhere. In their handy dispenser, CLINET® 101 are as easy to use as they are effective : CLINET® 101 wipes guarantee cleanness and cleanliness in a stroke.

EVERYDAY PROTECTION

Good hygiene is a daily battle. CLINET® 101 helps you to protect your environment and yourself. With CLINET® 101 on your side, you know that germs and dust just haven't got a chance !



interscience

The finishing touch

CLINET® 101

Laboratory





Cleaning and disinfection of laboratory equipment, glassware, benchtops, small instruments, keyboards, hands...

Hospital

Protection of patients, cleaning and disinfection of skin before injection, ambulances, waiting rooms, hands...

Office

Cleaning and disinfection of telephones, computer keyboards, washrooms, door handles, reception rooms, company cars...



Workshop

Cleaning and disinfection of instruments, tools, bench-tops, telephones, toilets, all frequently handled objects...

At Home

Protection of the family, cleaning and disinfection of taps, toilets, telephones, hands, door handles, car, thermometers...

COMPOSITION

Alcoholic solution of chlorhexidine digluconate (ethanol 50% p/p)
Bactericide : NF EN 1040 - NF T 72.151 - NF T 72.170 - NF T 72.190 -
Fungicide : NF EN 1275 - NF T 72.201 - NF T 72.190 - **Virucide** active on Herpes, Hepatitis B and HIV 1 (in vitro).

PRESENTATION

Handy dispenser with 120 pre-cut wipes (size : 13 x 20 cm). Ready to use : simply pull, tear and wipe. Boxes of 6 dispensers.
Reference : CLINET® 101.

INSTRUCTIONS FOR USE

Pull and tear off the wipe, clean the surface to be treated, allow to work without rinsing.

PRECAUTIONS

CLINET® 101 are hypoallergenic. For external use only - do not apply on mucous membranes, eyes, baby skin.

CLINET® 101 for everyday protection

Quand il faut que
ce soit propre



Vraiment propre



Please send me : a price offer for Clinet® 101.
 technical information on Clinet® 101.
 a free sample of Clinet® 101.

Company

Name

Position

Address

Zip code/City

Country

Tel. Fax

E-mail

Hygiene & Security

	info.	samples
Clinet® 202 (chlorhexidine individual cleaning wipes)	<input type="checkbox"/>	<input type="checkbox"/>
Clinet® 303 (bactericidal spray)	<input type="checkbox"/>	<input type="checkbox"/>
Clinet® anti-algae (algae-inhibiting liquid)	<input type="checkbox"/>	<input type="checkbox"/>
Anabac® (autoclave deodorant)	<input type="checkbox"/>	<input type="checkbox"/>
BioParfum® (odor neutralizers)	<input type="checkbox"/>	<input type="checkbox"/>

BagSystem® for sample preparation

BagMixer® (sample blenders)	<input type="checkbox"/>	
BagFilter® (sterile blender bags)	<input type="checkbox"/>	<input type="checkbox"/>
BagPipet® (automatic pipette)	<input type="checkbox"/>	

Please fax this request to your local dealer or to :

interscience

30 ch. Bois des Arpents 78860 St Nom - France
Tel. +33 (0)1 34 62 62 61 Fax +33 (0)1 34 62 43 03
E-mail : info@interscience.fr Web : www.interscience.fr