

Clinet® 202

Individually-packed cleaning wipes

interscience

Technical notice n° 2161 A (GB)

THE FINISHING TOUCH

Individually-packed Clinet® 202 cleaning wipes can be used wherever a high level of hygiene is needed. Their compact format allows them to be carried in a pocket, bag, car glove compartment... and are immediately available in any emergency.

CHARACTERISTICS

The disinfectant used is an alcoholic solution of benzalchonium chloride and chlorhexidine digluconate.

Renowned for its remarkable effectiveness, Clinet® 202 solution is hypoallergenic and complies to the following norms :

Bactericidal : NF EN 1040 -

NFT T 72.151 - NFT T 72.170 - NFT T 72.190

Fungicidal : NF EN 1275 - NF T 72.201 - NF T 72.190

Virucidal : Active on Herpes, Hepatitis B and HIV 1 (in vitro).

Clinet® 202 are ready for use. Simply wipe down the surface to be disinfected.

No need to rinse.

Individually-packed wipes.

Wipe size : 20 x 30 cm.


Box of 100 wipes.

1.4 kg - 21.5 x 26 x 10.5 cm.

Precautions :

For external use only.

Do not apply on mucous membranes or baby skin.



For cleaning and disinfecting large surfaces in the laboratory, production shop and all places of public affluence.

Production line, laboratory equipment, lavatory, telephone, ambulance, car, stretcher, public places, hands...

Please send me : a price offer for Clinet® 202.
 technical information on Clinet® 202.
 a free sample of Clinet® 202.

Company

Name

Position

Address

Zip code/City

Country

Tel. Fax

E-mail

Hygiene & Security

Clinet® 101 (chlorhexidine cleaning wipes) info. samples

Clinet® 303 (bactericidal spray)

Clinet® anti-algae (algae inhibiting liquid)

Anabac® (autoclave deodorant)

BioParfum® (odor neutralizers)

BagSystem® for sample preparation

BagMixer® (sample blenders)

BagFilter® (sterile blender bags)

BagPipet® (automatic pipette)

Please fax this request to your local dealer or to :

interscience

30 ch. Bois des Arpents 78860 St Nom - France

Tel. +33 (0)1 34 62 62 61 Fax +33 (0)1 34 62 43 03

E-mail : info@interscience.fr Web : www.interscience.fr