

Clinet® 303

Bactericidal spray

interscience

- Ready to spray
- Works on all surfaces, however uneven !
- No need to rinse



No contact
with hands
=
no risk of
irritation



Disinfects corners and
inaccessible areas

ALWAYS READY TO DISINFECT !

Clinet® 303 Spray is always near at hand. The active Clinet® 303 solution ensures highly effective disinfection and dries quickly leaving no trace. As a preventive measure at the end of the day, or in emergency cases, Clinet® 303 covers very large surfaces and reaches places inaccessible by wiping. Clinet® 303 guarantees excellent laboratory hygiene.

Bactericidal spray

TECHNICAL CHARACTERISTICS

Bactericidal. Fungicidal. AFNOR NFT 72.190.*

The Clinet® 303 solution is an alcohol-based preparation of glutaraldehyde, glyoxal and dimethyl-didecyl ammonium chloride. Clinet® 303 spray is always near at hand and is used PURE, by direct spraying on surfaces.

Do not rinse.

Net volume : 1 liter.

Dosage : ± 40 g/m² (well dampened).

Flash point < 22.8°C : inflammable.

Reference : Clinet® 303.

Box of 10 x 1 liter and 2 triggers.

*Technical file available on request.

Precautions : do not bring into contact with foodstuffs without rinsing (water); do not use on skin ; do not spray on methyl methacrylate surfaces (e.g. Perspex®), or if coated with alcohol-based varnish or oil-based paint ; keep away from all sources of heat.



CLINET® disinfectant line

Clinet® 101
impregnated chlorhexidine
disinfecting wipes

Clinet® 303
disinfecting spray

Clinet® 202 disinfecting wipes in individual sachets



Please send me : a price offer for Clinet® 303.
 technical information on Clinet® 303.
 a free sample of Clinet® 303.

Company
 Name
 Position
 Address

 Zip code/City
 Country
 Tel. Fax
 E-mail

Hygiene & Security

	info.	samples
Clinet® 101 (chlorhexidine cleaning wipes)	<input type="checkbox"/>	<input type="checkbox"/>
Clinet® 202 (chlorhexidine individual cleaning wipes)	<input type="checkbox"/>	<input type="checkbox"/>
Clinet® anti-algae (algae-inhibiting liquid)	<input type="checkbox"/>	<input type="checkbox"/>
Anabac® (autoclave deodorant)	<input type="checkbox"/>	<input type="checkbox"/>
BioParfum® (odor neutralizers)	<input type="checkbox"/>	<input type="checkbox"/>

BagSystem® for sample preparation

BagMixer® (sample blenders)	<input type="checkbox"/>	
BagFilter® (sterile blender bags)	<input type="checkbox"/>	<input type="checkbox"/>
BagPipet® (automatic pipette)	<input type="checkbox"/>	

Distributor

Please fax this request to your local dealer or to :